


FILTERS




5096HEPA



5096



303



303HEPA

- A family of multifunction filters available from one source.
- Each filter is designed to provide the best possible combination of:
 - ◆ High efficiency
 - ◆ Low compressible volume
 - ◆ Low resistance
 - ◆ Lightweight
 - ◆ Ergonomic
- Conical fittings meet ISO standard to ensure a secure connection.

Product Description	Bacterial/Viral filters for Anesthesia and Respiratory applications				
Ordering Info (50/case)	303 Respigard II™	303HEPA	5096	5096HEPA	30303*
Application	Machine-End Anes/Resp	Machine-End Anes/Resp	Patient-End Anes/Resp	Patient-End Anes/Resp	Bear 1™ & 2™, Resp
Filtration Efficiency					
BFE ¹	99.997%	99.9999%	99.997%	99.9999%	99.997%
VFE ¹	99.94%	99.999%	99.94%	99.999%	99.94%
Compressible Volume	30 mL	30 mL	27 mL	27 mL	30 mL
Minimum Tidal Volume	90 mL	90 mL	90 mL	90 mL	
Resistance to Flow @ 60 LPM	1.8 cm H ₂ O	2.2 cm H ₂ O	1.8 cm H ₂ O	2.2 cm H ₂ O	<2 cm H ₂ O
Connectors	22 mm F x 22 mm M / 15 mm F	22 mm F x 22 mm M / 15 mm F	22 mm F / 15 mm M x 22 mm M / 15 mm F	22 mm F / 15 mm M x 22 mm M / 15 mm F	22 mm F x 22 mm M / 15 mm F
Weight	16 g	16 g	17.3 g	18.5 g	16 g
Ports, Attachments	None	None	Gas Sampling Port	Gas Sampling Port	Mainline Kit

*Not pictured above.

Caution: Federal law restricts this device to sale by or on the order of a physician.

¹ Nelson Laboratories, Utah, USA. Data on file.

TD048 Rev.3 021104 © 2004 Vital Signs, Inc. All Rights Reserved.



20 Campus Road, Totowa, NJ 07512, USA • TEL: 800-932-0760 • FAX: +1-973-790-3307 • www.vital-signs.com
 Authorized European Representative: 13 -14 Eldon Way, Lineside Industrial Estate, Littlehampton, West Sussex, UK

