

Vital Signs
A GE Healthcare Company

CPAP Systems and
Accessories

**Comfort
Accuracy
and
High Flows**



imagination at work

Early Intervention Helps Decreases ICU Time



- **Flow Generator**

Variable Oxygen~ FiO₂ settings range from 33% - 100%
 High Flow Rates ~ Generates greater than 100 L/Min of flow to help reduce the work of breathing in acutely ill patients and potentially avoid intubation.

- **CPAP Mask**

Air cushion mask molds comfortably to the face; while providing a firm seal. Two sizes fit a broad range of patients, each with an integrated manometer port to verify PEEP settings and anti-asphyxia port.

- **AccuPEEP™ PEEP Valves**

Flow Independent ~ Threshold resistor valves are available in a wide range from 5 cmH₂O to 20 cmH₂O. The PEEP valves hold pressure constant, independent of flow.

Safety ~ Each PEEP valve comes with safety holes to help protect the patient in the event of exhalation valve occlusion.

Additional configurations are available. For more information please contact your local account representative.

Vital Signs, Inc.
 A GE Healthcare Company
 20 Campus Road
 Totowa, NJ 07512
 USA

www.vital-signs.com



Part Number	Description of Product kit / unit	U/M
C7111X31	CPAP Set 1.8 m hose, 303 filter, adult mask, 5 cmH ₂ O PEEP valve, silicone head strap, 72" (1.8 m) bubble tubing	10
C631XX12	CPAP Set 72" (1.8 m) hose, 303 filter, adult mask, 10 cmH ₂ O PEEP valve, silicone head strap	10
C6X1XX1X	CPAP Set 72" (1.8 m) hose, 303 filter, adult mask	10
9002C5	Adult CPAP Mask, Inlet 22 mm Male / Outlet 30 mm Female Inlet relief valve, 9260 head strap	10
9002	Adult CPAP Mask, Inlet 31 mm Female / Outlet 30 mm Female, One-way valves on inlet and outlet ports	10
9003B	2.5 cmH ₂ O fixed Accu-PEEP valve	10
9005B	5.0 cmH ₂ O fixed Accu-PEEP valve	10
9007B	7.5 cmH ₂ O fixed Accu-PEEP valve	10
9010B	10.0 cmH ₂ O fixed Accu-PEEP valve	10
9012B	12.5 cmH ₂ O fixed Accu-PEEP valve	10
9015B	15.0 cmH ₂ O fixed Accu-PEEP valve	10
9022B	20.0 cmH ₂ O fixed Accu-PEEP valve	10
9260	Silicone Head strap	20
9100	Head strap, Silicone and Polyester	10

© 2009 General Electric Company - All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE, GE Monogram, AccuPEEP, and Vital Signs are trademarks of General Electric Company.

Vital Signs, Inc., a General Electric company, doing business as GE Healthcare.

L2007 CPAP Broch EN Rev. Rel. 05/2009

ACCC-0294-09.09-EN-US