

# PRE-HOSPITAL CPAP SYSTEMS



- Flows greater than 100 lpm
- O<sub>2</sub> concentrations from 33% - 100%
- Lightweight and portable
- Attaches to any 50 psi connector

Product Description	Portable flow & CPAP system for pre-hospital care of Congestive Heart Failure patients		
Ordering Information	Product Code	Description	Qty/Case
	9250	Downs Adjustable FIO <sub>2</sub> Flow Generator	1
	6098GAP	Airway Pressure Gauge with Gauge Protector	1
	C631XXX1	72" Circuit, CPAP Mask, 10 cm PEEP Valve, Headstrap	10
	C631XX11	72" Circuit, CPAP Mask, 10 cm PEEP Valve, Headstrap and Filter	10
	C731XXX1	72" Circuit, CPAP Mask, 10 cm PEEP Valve, Headstrap, Manometer Line	10
	C731XX11	72" Circuit, CPAP Mask, 10 cm PEEP Valve, Headstrap, Manometer Line and Filter	10
Expected Outcomes	Higher blood oxygen levels with lower FIO <sub>2</sub> s. Lower incidence of intubation upon arrival at hospital.		

20 Campus Road, Totowa, NJ 07612, USA • TEL: 800-932-0760 • FAX: +1-973-790-3307 • www.vital-signs.com  
 Sussex Business Village, Lake Lane, Barnham, PO22 0AL, UK • TEL: +44-1243-555300 • FAX: +44 1243-555400 • vitalcare@vital-signs.co.uk

**vital signs inc.**

**vital**

