

THE BERMAN INTUBATING/PHARYNGEAL AIRWAY



**LATEX
FREE**

- Ideal for fiberoptic bronchoscopy and gastroscopic procedures for guiding, positioning, and protecting costly fiberoptic equipment.
- Large internal diameter accommodates the passage of an endotracheal tube through its center.
- May be opened longitudinally and separated from the tracheal tube.

Product Description	Intubating/Pharyngeal airway		
Ordering Information (50/case)	Product Code	Description	Size
	2408	No. 8 Small (Child)	8 cm
	2409	No. 9 Medium (Adult)	9 cm
	2410	No. 10 Large (Adult)	10 cm

20 Campus Road, Totowa, NJ 07512, USA • TEL: 800-932-0760 • FAX: +1-973-790-3307 • www.vital-signs.com
 Authorized European Representative: Vital Signs, Ltd., 13-14 Eldon Way, Lineside Industrial Estate, Littlehampton, West Sussex, UK